

A CLEANER & DISINFECTANT PRODUCT FOR EVAPORATOR COILS IN THE HVACR INDUSTRY

- Cleans & Deodorizes in one easy step
- EPA registered for the HVAC application
- Ready-to-use disinfectant formulation
- No rinsing required
- Eliminates odors
- Cuts dirt and grime
- Can be used with Talon (4776-0)

Description

Evap-Fresh No Rinse is a ready-to-use, non-abrasive, one-step cleaner and disinfectant for evaporator coils in small air conditioners, commercial air handling equipment, HVAC cooling coils and refrigeration equipment. It is the first product of its kind in the industry to receive EPA registration for such claims. For hard to reach areas, use in conjunction with 27 inch Talon Flexible Spray Wand.

Evap-Fresh has been tested effective via AOAC methodology against a wide variety of organisms:

- Salmonella typhi
- Yersinia enterocolitica
- Enterococcus faecium
- Corynebacterium ammoniagenes
- Methicillin resistant Staphylococcus aureus (MRSA)
- Vancomycin resistant Enterococcus faecalis
- Methicillin resistant Staphylococcus epidermidis
- Vancomycin intermediate resistant Staphylococcus aureus
- Pseudomonas aeruginosa
- Staphylococcus aureus
- Salmonella choleraesuis
- Escherichia coli
- Streptococcus pyogenes
- Listeria monocytogenes

Indoor Air Quality Products

Evap-Fresh™ No Rinse



Directions for Use

Gallon:


Remove heavy soil prior to product use. Pour product into a Nu-Calgon sprayer of choice and then hold sprayer six to eight inches from surface to be treated. Spray area until it is covered with the solution. Allow product to penetrate and remain wet for 3 minutes. No scrubbing necessary. The product will not leave grit or soap scum.

Aerosol:

Shake can before using. Hold container six to eight inches from surface to be treated. Spray area until it is covered with white foam. Allow product to penetrate and remain wet for 10 minutes. No scrubbing necessary. The product will not leave grit or soap scum. Lemon scented.

Packaging

- | | |
|----------------------|----------------|
| 1 gallon | 4166-08 |
| 18 oz. aerosol spray | 4166-75 |



Additional Information

- Exhibits effective disinfectant activity against *Clostridium difficile* vegetative form, AATC 9689, when the treated surface is allowed to remain wet for 1 minute (60 seconds).
- Tuberculocidal Activity - Evap-Fresh exhibits disinfectant efficacy against *Mycobacterium tuberculosis* (BCG) at 20 degrees centigrade when the treated surface is allowed to remain wet for 5 minutes.
- Virucidal Activity - Evap-Fresh kills on hard, non-porous inanimate surfaces when allowed to remain wet for a 10 minute contact time against: Hepatitis A Virus (HAV), Canine Parvovirus and Poliovirus Type 1.
- When used on environmental, inanimate hard surfaces exhibits virucidal activity against Hepatitis B Virus (HBV), Hepatitis C Virus (HCV), and Bovine Viral Diarrhea Virus (BVDV) when the surface is allowed to remain wet for 5 minutes.
- When used on environmental, inanimate hard surfaces exhibits virucidal activity against Human Coronavirus and SARS associated Coronavirus when the treated surface is allowed to remain wet for 2 minutes.
- When used on environmental, inanimate hard surfaces exhibits virucidal activity against Rhinovirus type 39 (ATCC VR-340) and Rotavirus when the treated surface is allowed to remain wet for 3 minutes.
- When used on environmental, inanimate hard surfaces exhibits virucidal activity against HIV-1 (associated with AIDS) when the treated surface is allowed to remain wet for 1 minute (60 seconds).
- When used on environmental, inanimate hard surfaces exhibits effective virucidal activity against Norovirus (Norwalk virus), Feline Calicivirus and Rabies Virus when the treated surface is allowed to remain wet for 30 seconds.
- Kills HIV-1, HBV and HCV on pre-cleaned surfaces/objects

Read and understand the product's label and Safety Data Sheet
("SDS") for precautionary and first aid information.

The SDS is available on the Nu-Calgon website at www.nucalgon.com.

